




Exploring the innovative world of dental materials and equipment



Website: https://jeedent.com/

E-mail: Eunice@jeedent.com/Sales@jeedent.com

Whatsapp:+86 18968015464











Company Introduction

Our Story

Jeedent is Hangzhou Jideng Medical
Technology Co., Ltd., is a professional
manufacturer of dental products integrating
R&D, production, and sales, providing support
to dental laboratories and clinics worldwide.
We have a professional sales and service team
to better serve our global customers with
expert guidance and tailored solutions.







Factory Introduction

Producing Department



Regulated Production



Printing



Packaging



Shipping





Dental Consumables Classification

- 1. Stain & Glaze
- 2. Dental Milling Burs
- 3. Dental polish kit
- 4. Sintering Paste
- **5. Polishing Paste**
- 6. Scanning Gel / Scan Powder





Stain & Glaze (1)



Set Composition:

Includes 4 primary ABCD shades, 18 special effect colors, glaze, 3D natural shade, and diluting & refreshing liquid

Scope of application:

Full Zircon, Copings, PFM, Glass ceramic

Stain & Glaze (2)



Core Advantage:

- 1. Each paste bottle contains glaze ingredients, no separate glaze application required
- 2. Simplified operation only one stain application needed to achieve the desired effect
- 3. Glaze paste contains fine particles, no powder-liquid mixing required before use
- 4. Flexible crown shading paste shades can be mixed as desired, ready for direct use after even blending



Dental Milling Burs (1)

Application Scenarios:

- 1. Dental Restoration Fabrication
- 2. Tooth Preparation
- 3. Oral Surgery
- 4. Dental Material Processing





Dental Milling Burs (2)



Core Advantages:

- 1. High-Precision Machining
- 2. Exceptional Durability
- 3. Efficient Cutting
- 4. Reduced Thermal Damage
- 5. Strong Compatibility
- 6. Superior Surface Finish



Dental polish kit (1)



Key Uses

- **1.** Restore Polishing: Smooth ceramic crowns, composite fillings, or resin veneers post-fabrication/bonding.
- 2. Natural Teeth Care: Remove stains (coffee/tea) and polish enamel (ideal after whitening).
- **3.** Orthodontic Aftercare: Smooth enamel post-brace removal; relieves discomfort.
- **4.** Lab Finishing: Final touch for inlays/onlays/denture bases to meet standards.



Dental polish kit (2)



Core Advantages

- **1.** All-in-One Solution--Includes tools (pastes, discs, brushes) for various dental materials (ceramics, composites, metals) no extra purchases needed.
- **2.** Pro-Grade Finish--High-purity abrasives deliver smooth, glossy surfaces; reduces bacteria adhesion and boosts aesthetics.
- **3. Safe & Compatible**--Matches material hardness to avoid scratching; meets biocompatibility standards (no oral irritation).
- 4. Easy & Efficient -- Ergonomic parts + clear guidelines save time for dentists/lab techs.
- **5.** Long-Lasting--Durable components cut replacement costs.



Sintering Paste



Basic Information

1. Specifications and Models: 7g, 14g, 42g, 70g

2. Color: White

3. Form: In room temperature conditions, the paste form

Advantages

1. Easy to remove, no residue

2. Provides strong support, ensuring stable operation





Applicable Restorations:

Veneers, inlays, onlays, partial crowns, single crowns, pontics

Core Function 1:

Stabilizes & supports restorations during porcelain firing or crystallization

Core Function 2:

Prevents restorations from deforming in the crystallization process



Polishing Paste



Specification	10, 30g
Mode1	ZS25, SC25, ZS50, SC50
XX can be customized	ZSXX, SCXX

Product Application Scope

Used outside the mouth to grind and polish the restoration, making its surface smooth and uniform

Advantages

- 1. Efficient polishing, significantly enhancing brightness
- 2. No glazing, simplifying operation procedures
- 3. excellent durability performance



Scanning Gel / Scan Powder



Specification	10, 20g
Color	White
Form	Creamy Paste, Powder

Product Application Scope

High-precision scanning of precise components such as implants

Advantages

- 1. Easy to clean, splash-proof, durable
- 2. Effectively prevents dust inhalation, protects health
- 3. transportation safety, easy to operate







Shipping

To better ensure the safety of your goods, we will provide professional, environmentally friendly and efficient packaging services. We directly cooperate with DHL, UPS and FedEx, and send your products to you with fast and secure services.












Certificate

We comply with CE 0197 and FDA regulations and are certified to ISO 9001 and ISO 13485 quality management standards.





Registration No.

Hangzhou Jideng Medical Technology Co., Ltd. Manufacturer:

Room 601, 6th Floor, Building 1, No. 96 Fangxingdu Street,

Xingqiao Subdistrict, Linping District,

Hangzhou City,

Zhejiang CN-MF-000025755

EUDAMED Single Registration No.:

Products:

Products of class IIa: Q010101 - DENTAL RESTORATION DEVICES

Zirconia Dental Ceramics

Authorized representative(s): Riomavix S I

Calle deAlmansa 55.1D.Madrid 28039 Spain

Annex XI Part A

Revision: Description Issue date 2025-07-23 Initial certification

The Notified Body bereby declares that the requirements of Annex XI, Part A of the REGULATION (EU) 2017/745 have been met for the listed products. The above named manufacturer has established and applies a production quality assurance, which is subject to periodic surveillance, defined by Annex XI, Part A, Section 7 of the aforementioned regulation.

If class III devices or class IIb devices are covered by this certificate, an EU type-examination certificate in accordance

with Annex X of the aforementioned regulation is required before placing them on the marke

244451789-200

2025-07-23 Effective date: 2030-07-22 Issue date: 2025-07-23

> TÜV Rheinland LGA Products GmbH Tillystraße 2 · 90431 Nürnberg · Germany

TÜV Rheinland LGA Products GmbH is a Notified Body according to REGULATION (EU) 2017/745 concerning medical









Thank you for your attention!

Website: https://jeedent.com/

E-mail: Eunice@jeedent.com/Sales@jeedent.com

Whatsapp:+86 18968015464







